


OSCAR VISION

From Darkness to Light



VISUAL ACUITY

FINE VISION CHART - AB6

Advance Features Fine Vision Chart - AB6

LIST OF CHARTS

- Landolt "C" Ring
- Maddox Test
- Tumbling E
- Peripheral Aquity
- Alphabets (Random)
- Red Green (or) Duochrome
- Numeric (Random)
- Cross Cylinder Test
- Dots Chart
- Ishihara Colour Charts
- Pediatric Charts
- Animal Parade Charts
- Contrast Sensitivity
- Cartoon Charts
- Logmar Charts
- Astig Fan
- Snellen Charts
- All Reginal Language
- ETDRS Contrast Charts, Urdu & Arabic

Chart Distance Available :

6 feet, 8 feet, 10 feet, 12 feet, 15 feet, 20 feet and Mirror Charts. All distance charts in built with machine

Monitor Size : 18.5 inch (Wide)

Eye Educational Chart:

Cataract, Retina, Cornea, Glaucoma, Lasik, Contact Lens (Total 100 Charts)

Voltage : 110/230V/50 or 60 HZ ; Weight : 3 kg Net

All Specification subject to change without notice

Your Nearest Dealer :

Contact us :

OSCAR VISION TECH

71/2, 12th Avenue, R.R. Colony, 1st Street,
Ahok Nagar, Chennai - 600 083, INDIA

Phone : +91 44 23710690

Mobile : +91 94444 13257

Email : oscarjayakumar@yahoo.in
oscarjayakumar@gmail.com