

OPMI LUMERA 300 from ZEISS



OPMI LUMERA 300 from ZEISS



Technical Data and Configurations

OPMI LUMERA 300 from ZEISS

Configure the ZEISS OPMI LUMERA 300 surgical microscope to your individual needs. Select from a variety of accessories such as an integrated HD camera or an assistants package with an assistant's microscope that doesn't take any light from the surgeons microscope. Whatever you select you can always be confident to get best-in-class visualization from ZEISS.

	Basic	Standard	Advanced
Optic			
Focus range	48 mm	48 mm	48 mm
Objective lens	f = 200 mm	f = 200 mm	f = 200 mm
Wide-angle eyepiece	10x	12.5x	12.5x
Magnification	Manual 5-step magnification changer	Motorized zoom	Motorized zoom
Binocular tube	45° Inclined tube f = 170 mm	45° Inclined tube f = 170 mm	180° Swivel tube f = 170 mm

OPMI LUMERA 300 from ZEISS

Apochromatic optic



Integrated assistant interface

